MAGERI benchmark using reference standard DNA library

Mikhail Shugay

December 1, 2016

```
Load metadata
df.vmeta <- read.table("hd734_variant_metadata.txt", sep="\t", header=T) # variants observed in HD734 a
df.smeta <- read.table("sample_metadata.txt", sep="\t", header=T) # metadata for amplicon sequencing sa
VCF parsing function
library(stringr)
library(pROC)
## Type 'citation("pROC")' for a citation.
## Attaching package: 'pROC'
## The following objects are masked from 'package:stats':
##
##
       cov, smooth, var
library(ggplot2)
library(ggbeeswarm)
library(dplyr)
##
## Attaching package: 'dplyr'
## The following objects are masked from 'package:stats':
##
##
       filter, lag
## The following objects are masked from 'package:base':
##
##
       intersect, setdiff, setequal, union
read_vcf <- function(file_name) {</pre>
  .vcf <- read.table(file_name, header = F, sep = "\t", stringsAsFactors = F)</pre>
  colnames(.vcf) <- c("chromosome", "position", "skip1", "from", "to", "qual", "skip2", "info", "skip3"</pre>
  .vcf$skip1 <- NULL</pre>
  .vcf$skip2 <- NULL
  .vcf$skip3 <- NULL
  .vcf$skip4 <- NULL
  .vcf$qual <- as.integer(.vcf$qual)</pre>
  .vcf <- subset(.vcf, nchar(from) == 1 &</pre>
                    nchar(to) == 1 & !is.na(qual)) # no indels
  .infosplit <- str_split_fixed(.vcf$info, regex("[=;]"), 15)[,c(2, 4, 12, 14)]</pre>
  .vcf$coverage <- as.numeric(.infosplit[,1])</pre>
```

```
.vcf$frequency <- as.numeric(.infosplit[,2])</pre>
  .vcf$bb.a <- as.numeric(.infosplit[,3])</pre>
  .vcf$bb.b <- as.numeric(.infosplit[,4])</pre>
  .vcf$info <- NULL</pre>
  .vcf$count <- as.integer(round(.vcf$coverage * .vcf$frequency))</pre>
  .vcf
}
head(read_vcf("p126.h4_2_ballast_m1.vcf"))
## Warning in read_vcf("p126.h4_2_ballast_m1.vcf"): NAs introduced by coercion
##
     chromosome position from to qual coverage
                                                     frequency
                                            2265 0.0004415011 0.9957058
## 1
           chr2 212295704
                              C A
                                     8
## 2
           chr2 212295705
                              C A
                                      8
                                             2265 0.0004415011 0.9957058
## 3
           chr2 212295713
                              A G
                                      7
                                            2265 0.0004415011 1.0295689
## 4
           chr2 212295718
                            G A 6
                                         2265 0.0004415011 2.1314828
## 5
           chr2 212295725
                           C A 8
                                            2264 0.0004416961 0.9957058
## 6
           chr2 212295732 C T
                                     24
                                             2264 0.0013250883 2.1314828
##
         bb.b count
## 1 21862.23
## 2 21862.23
## 3 16394.12
                  1
## 4 24344.82
                  1
## 5 21862.23
                  1
## 6 24344.82
Read samples with HD734 standard DNA and control human DNA, append metadata
library(TailRank) # For betabinom.
## Loading required package: oompaBase
# installing:
# source("https://bioconductor.org/biocLite.R")
# biocLite("Biobase")
# install.packages("TailRank", repos="http://R-Forge.R-project.org")
df <- data.frame()</pre>
read_vcf_with_metadata <- function(file_name, primer_set, replica, ratio, type) {</pre>
  .vcf <- read_vcf(file_name)</pre>
  .vcf <- merge(.vcf, df.vmeta, all.x = type != "standard", all.y = F)</pre>
  .vcf$known.frequency <- .vcf$known.frequency * ratio</pre>
  .vcf\$known.frequency[is.na(.vcf\$known.frequency)] <- 0
  .vcf$primer_set <- primer_set</pre>
  .vcf$replica <- primer_set</pre>
  .vcf$type <- type</pre>
  .vcf <- subset(.vcf, frequency < 0.4 & count > 0) # remove alleles in control
  .vcf
}
```

```
for (i in 1:nrow(df.smeta)) {
  df <- with(df.smeta, rbind(df,</pre>
                             read_vcf_with_metadata(paste(prefix[i], "vcf", sep="."),
                                                    primer set[i],
                                                    replica[i],
                                                    ratio[i],
                                                    type[i])))
}
## Warning in read_vcf(file_name): NAs introduced by coercion
# Compute scores
#df$qual <- with(df,
\# mapply(function(x,y,a,b) \ 1 - pbb(x,y,a,b) + 0.5 * dbb(x,y,a,b),
         count, coverage, bb.a, bb.b)
# )
```

```
#df$qual <- ifelse(df$qual <= 1e-100, 999, -10*log10(df$qual))
Group observed variants into tiers
df$tier <- cut(df$known.frequency, c(-1, 0, 0.009, 0.02, 1))
levels(df$tier) <- c("error", "0.1%", "1%", "5%")</pre>
summary(df$tier)
## error 0.1%
                  1%
                         5%
    2468
##
            43
                  46
                         12
Number of detected variants and sampling
df.v.summary <- df.vmeta</pre>
df.v.summary$known.frequency <- df.v.summary$known.frequency * 0.1 # 1:9 diluted samples
df.v.summary$ratio <- "1:10"
df.v.summary.2 <- df.vmeta
df.v.summary.2$ratio <- "1:1"</pre>
df.v.summary <- rbind(df.v.summary, df.v.summary.2)</pre>
df.v.summary$tier <- cut(df.v.summary$known.frequency, c(-1, 0, 0.009, 0.02, 1))
levels(df.v.summary$tier) <- c("error", "0.1%", "1%", "5%")</pre>
print(df.v.summary)
##
                                                                   id
      chromosome position from to
## 1
            chr7 140453136
                                               BRAF V600E, BRAF V600K
## 2
            chr7 140453136
                               A C
                                               BRAF_V600G, BRAF_V600R
## 3
                                  T BRAF_V600M, BRAF_V600K, BRAF_V600R
            chr7 140453137
                               C
## 4
            chr7 55241707
                              G A
                                                           EGFR G719S
## 5
            chr7 55259515
                              T G
                                                           EGFR L858R
            chr7 55259524
## 6
                               T A
                                                           EGFR L861Q
## 7
            chr7 55249071
                               С
                                 Τ
                                                           EGFR T790M
                               A T
## 8
           chr19
                   3118942
                                                          GNA11_Q209L
## 9
            chr2 209113112
                               C T
                                                           IDH1_R132C
## 10
            chr2 209113113
                               G A
                                                           IDH1_R132H
## 11
           chr12 25378562
                               C
                                  Τ
                                                           KRAS_A146T
## 12
           chr12 25398284
                                 G
                                                            KRAS_G12A
## 13
                  25398285
                               C A
                                                            KRAS_G12C
           chr12
## 14
           chr12
                  25398284
                               С
                                  Τ
                                                            KRAS_G12D
                                 G
## 15
           chr12 25398285
                               С
                                                            KRAS_G12R
## 16
           chr12 25398285
                               C T
                                                            KRAS G12S
           chr12 25398284
                               C A
                                                            KRAS_G12V
## 17
## 18
           chr12
                  25398281
                               С
                                  Т
                                                            KRAS_G13D
## 19
           chr12 25380275
                               Т
                                  G
                                                            KRAS_Q61H
                                                            KRAS_Q61L
## 20
           chr12 25380276
                               Τ
                                 Α
            chr1 115256528
## 21
                               Т
                                                            NRAS Q61H
                                 Α
## 22
            chr1 115256530
                               G
                                  Τ
                                                            NRAS Q61K
## 23
            chr1 115256529
                               T A
                                                            NRAS_Q61L
## 24
            chr1 115256529
                               T C
                                                            NRAS_Q61R
## 25
            chr3 178936082
                               G
                                  Α
                                                         PIK3CA_E542K
## 26
            chr3 178936091
                               G
                                  Α
                                                         PIK3CA_E545K
## 27
                               A G
            chr3 178952085
                                                        PIK3CA_H1047R
## 28
            chr7 140453136
                               Α
                                  Τ
                                               BRAF_V600E, BRAF_V600K
```

T BRAF_V600M, BRAF_V600K, BRAF_V600R

BRAF_V600G, BRAF_V600R

EGFR_G719S

29

30

31

chr7 140453136

chr7 140453137

chr7 55241707


A C

С

```
## 32
            chr7 55259515
                                                          EGFR L858R
## 33
                              T A
            chr7 55259524
                                                          EGFR_L861Q
                              C T
## 34
            chr7
                  55249071
                                                          EGFR T790M
                                                         GNA11_Q209L
## 35
           chr19
                   3118942
                                 Τ
                              Α
## 36
            chr2 209113112
                                 Τ
                                                          IDH1 R132C
## 37
            chr2 209113113
                              G A
                                                          IDH1 R132H
                                                          KRAS_A146T
## 38
                  25378562
                                 Т
           chr12
                              С
                                                           KRAS_G12A
## 39
           chr12
                  25398284
                                 G
## 40
           chr12
                  25398285
                              C A
                                                           KRAS G12C
## 41
                              C T
           chr12
                  25398284
                                                           KRAS_G12D
## 42
           chr12
                  25398285
                              C G
                                                           KRAS_G12R
                              С
                                 Т
## 43
                  25398285
                                                           KRAS_G12S
           chr12
## 44
           chr12 25398284
                              C
                                 Α
                                                           KRAS_G12V
## 45
                  25398281
                                 Τ
           chr12
                                                           KRAS_G13D
## 46
                  25380275
                              T
                                 G
                                                           KRAS_Q61H
           chr12
## 47
           chr12
                  25380276
                              T
                                 Α
                                                           KRAS_Q61L
## 48
            chr1 115256528
                              T
                                 Α
                                                           NRAS_Q61H
## 49
            chr1 115256530
                                                           NRAS Q61K
## 50
            chr1 115256529
                              T A
                                                           NRAS_Q61L
## 51
            chr1 115256529
                              T C
                                                           NRAS Q61R
                              G A
## 52
            chr3 178936082
                                                        PIK3CA_E542K
## 53
            chr3 178936091
                              G A
                                                        PIK3CA E545K
## 54
            chr3 178952085
                              A G
                                                       PIK3CA_H1047R
      known.frequency reference ratio tier
##
## 1
                        BRAF E15 1:10 0.1%
               0.0080
               0.0010
                        BRAF E15 1:10 0.1%
## 3
               0.0010
                        BRAF_E15
                                 1:10 0.1%
## 4
               0.0167
                        EGFR_E18 1:10
                                          1%
## 5
               0.0010
                        EGFR_E21 1:10 0.1%
## 6
               0.0010
                        EGFR_E21 1:10 0.1%
## 7
               0.0010
                        EGFR_E20
                                  1:10 0.1%
## 8
               0.0013
                        GNA11_E5 1:10 0.1%
## 9
               0.0013
                         IDH1_E4 1:10 0.1%
## 10
               0.0013
                         IDH1_E4 1:10 0.1%
## 11
               0.0013
                         KRAS E4 1:10 0.1%
## 12
               0.0013
                         KRAS_E2 1:10 0.1%
## 13
               0.0013
                         KRAS E2 1:10 0.1%
## 14
               0.0013
                         KRAS_E2 1:10 0.1%
## 15
               0.0013
                         KRAS_E2 1:10 0.1%
## 16
               0.0013
                         KRAS_E2 1:10 0.1%
## 17
               0.0013
                         KRAS E2 1:10 0.1%
## 18
               0.0250
                         KRAS E2 1:10
                                          5%
               0.0013
                         KRAS E3 1:10 0.1%
## 19
## 20
               0.0013
                         KRAS_E3 1:10 0.1%
## 21
               0.0013
                         NRAS_E3 1:10 0.1%
## 22
                         NRAS_E3 1:10 0.1%
               0.0013
## 23
               0.0013
                         NRAS_E3
                                  1:10 0.1%
## 24
               0.0013
                         NRAS_E3
                                  1:10 0.1%
## 25
               0.0013 PIK3CA_E9
                                  1:10 0.1%
## 26
               0.0013
                       PIK3CA_E9
                                  1:10 0.1%
## 27
               0.0300 PIK3CA_E20
                                  1:10
                                          5%
## 28
               0.0800
                        BRAF_E15
                                   1:1
                                          5%
## 29
               0.0100
                        BRAF E15
                                   1:1
                                          1%
## 30
               0.0100
                        BRAF E15
                                   1:1
                                          1%
```

```
## 31
               0.1670
                        EGFR E18
                                          5%
                                    1:1
                        EGFR_E21
## 32
               0.0100
                                          1%
                                    1:1
## 33
               0.0100
                        EGFR E21
                                    1:1
                                          1%
## 34
               0.0100
                        EGFR_E20
                                    1:1
                                          1%
## 35
               0.0130
                        GNA11 E5
                                    1:1
                                          1%
## 36
               0.0130
                         IDH1 E4
                                   1:1
                                          1%
## 37
                         IDH1 E4
               0.0130
                                   1:1
                                          1%
                         KRAS_E4
## 38
               0.0130
                                    1:1
                                          1%
## 39
               0.0130
                         KRAS E2
                                    1:1
                                          1%
                                          1%
## 40
               0.0130
                         KRAS_E2
                                    1:1
## 41
               0.0130
                         KRAS_E2
                                    1:1
                                          1%
                         KRAS_E2
## 42
                                          1%
               0.0130
                                    1:1
## 43
               0.0130
                         KRAS_E2
                                    1:1
                                          1%
## 44
               0.0130
                         KRAS_E2
                                    1:1
                                          1%
## 45
               0.2500
                         KRAS_E2
                                    1:1
                                          5%
## 46
               0.0130
                         KRAS_E3
                                    1:1
                                          1%
## 47
                         KRAS_E3
                                          1%
               0.0130
                                    1:1
## 48
               0.0130
                         NRAS E3
                                    1:1
                                          1%
## 49
               0.0130
                         NRAS_E3
                                    1:1
                                          1%
## 50
               0.0130
                         NRAS E3
                                    1:1
                                          1%
## 51
               0.0130
                         NRAS_E3
                                    1:1
                                          1%
## 52
               0.0130 PIK3CA E9
                                          1%
                                    1:1
## 53
               0.0130 PIK3CA_E9
                                          1%
                                    1:1
## 54
               0.3000 PIK3CA E20
                                    1:1
                                          5%
df.v.summary <- df.v.summary %>%
  group_by(tier) %>%
  summarise(total = n())
df.v.summary$total <- 2 * df.v.summary$total # we have 2 independent experiments
df.v.summary <- merge(df.v.summary,</pre>
                      df %>% group_by(tier) %>%
                         summarise(found = n(),
                                   detect.prob.m = mean(1-exp(-frequency * coverage)),
                                   detect.prob.sd = sd(1-exp(-frequency * coverage))))
df.v.summary$detection.prob.lb <- with(df.v.summary,</pre>
                                        detect.prob.m - 1.96*detect.prob.sd - 5e-3)
df.v.summary$detection.prob.ub <- with(df.v.summary,</pre>
                                        detect.prob.m + 1.96*detect.prob.sd)
df.v.summary$detection.prob.ub <- ifelse(df.v.summary$detection.prob.ub > 1,
                                          1,
                                          df.v.summary$detection.prob.ub )
library(scales)
ggplot(df.v.summary, aes(x=tier)) +
  geom_bar(stat="identity", aes(y=detection.prob.ub), fill = "grey70", width=0.5)+
  geom_bar(stat="identity", aes(y=detection.prob.lb), fill = "white")+
  geom_point(aes(y = found/total), size = 5, color="red", shape="+") +
  geom_text(aes(y = found/total - 0.05, label = paste(found, total, sep="/"))) +
  scale_y_continuous("Detection rate", labels=percent) +
  xlab("Variant tier") +
  theme_bw() +
```






Quality score and error P-values

```
df.e <- df %>%
  filter(type == "blank") %>%
  mutate(true.p.value = -10 * log10(1 - (rank(frequency) - 0.5) / n()))


ggplot(df.e, aes(true.p.value, qual)) +
  geom_abline(slope = 1, intercept = 0, linetype = "dashed") +
  geom_smooth(color="red", fill="grey", alpha=0.5) +
  scale_x_continuous("-10 log10 P-value", limits = c(0, 40)) +
  scale_y_continuous("Q score", limits = c(0, 40)) +
  theme_bw()
```

Warning: Removed 1 rows containing non-finite values (stat_smooth).




Frequency of variants and errors

```
ggplot(df, aes(x=frequency, color = tier)) +
  geom_density(aes(fill=tier), alpha=0.7) +
  ylab("Density") +
  scale_x_log10("Frequency") +
  scale_color_brewer("Variant tier", palette = "Set1") +
  scale_fill_brewer("Variant tier", palette = "Set1") +
  theme_bw()
```



Variant quality scores

```
ggplot(df,aes(tier, qual, fill=tier)) +
  geom_hline(yintercept = 13.0103, linetype ="dotted") +
  geom_hline(yintercept = 20.0103, linetype ="dashed") +
  geom_quasirandom(varwidth = T, shape=21, color="grey10") +
  annotate("text", label = "p < 0.05", x=4.3,y=13.0103+2) +
  annotate("text", label = "p < 0.01", x=4.3,y=20.0103+4) +
  scale_fill_brewer("", palette = "Set1", guide=F) +
  xlab("Variant tier") + ylab("Q score") +
  theme_bw()</pre>
```



ROC curve for rare variants

```
df.1 <- subset(df, tier %in% c("error", "0.1%"))</pre>
df.1$type <- ifelse(df.1$tier == "error", 0, 1)</pre>
make_roc <- function(rocobj, type) {</pre>
  .df.roc \leftarrow data.frame(spec = seq(0, 1, 0.01))
  sens.ci <- ci.se(rocobj, specificities = .df.roc$spec)</pre>
  .df.roc$sens.lo <- sens.ci[,1]</pre>
  .df.roc$sens.me <- sens.ci[,2]</pre>
  .df.roc$sens.hi <- sens.ci[,3]</pre>
  .df.roc$type <- type</pre>
  .df.roc
}
rocobj <- roc(type ~ qual, df.1, ci=T)</pre>
print(rocobj)
##
## Call:
## roc.formula(formula = type ~ qual, data = df.1, ci = T)
## Data: qual in 2468 controls (type 0) < 43 cases (type 1).
## Area under the curve: 0.9277
## 95% CI: 0.8727-0.9827 (DeLong)
```

```
df.roc <- make_roc(rocobj, "Q score")</pre>
rocobj <- roc(type ~ frequency, df.1, ci=T)</pre>
print(rocobj)
##
## Call:
## roc.formula(formula = type ~ frequency, data = df.1, ci = T)
## Data: frequency in 2468 controls (type 0) < 43 cases (type 1).
## Area under the curve: 0.8637
## 95% CI: 0.7839-0.9435 (DeLong)
df.roc <- rbind(df.roc, make_roc(rocobj, "Frequency"))</pre>
ggplot(df.roc, aes(x=spec)) +
  geom_ribbon(aes(ymin=sens.lo, ymax=sens.hi, group=type), fill="grey", alpha=0.5) +
  geom_abline(slope = 1, intercept = 1, linetype = "dashed") +
  geom line(aes(y=sens.me, color = type)) +
  scale_x_reverse("Specificity") +
  scale_y_continuous("Sensitivity", limits=c(0,1)) +
  scale_color_brewer("Threshold", palette = "Set1") +
  theme_bw()
```

