MAGERI benchmark using reference standard DNA library

Mikhail Shugay March 1, 2017

```
Load metadata
df.vmeta <- read.table("hd734_variant_metadata.txt", sep="\t", header=T) # variants observed in HD734 a
df.smeta <- read.table("sample_metadata.txt", sep="\t", header=T) # metadata for amplicon sequencing sa
VCF parsing function
library(stringr)
library(pROC)
## Type 'citation("pROC")' for a citation.
## Attaching package: 'pROC'
## The following objects are masked from 'package:stats':
##
##
       cov, smooth, var
library(ggplot2)
library(ggbeeswarm)
library(dplyr)
##
## Attaching package: 'dplyr'
## The following objects are masked from 'package:stats':
##
##
       filter, lag
## The following objects are masked from 'package:base':
##
##
       intersect, setdiff, setequal, union
read_vcf <- function(file_name) {</pre>
  .vcf <- read.table(file_name, header = F, sep = "\t", stringsAsFactors = F)</pre>
  colnames(.vcf) <- c("chromosome", "position", "skip1", "from", "to", "qual", "skip2", "info", "skip3"</pre>
  .vcf$skip1 <- NULL</pre>
  .vcf$skip2 <- NULL
  .vcf$skip3 <- NULL
  .vcf$skip4 <- NULL
  .vcf$qual <- as.integer(.vcf$qual)</pre>
  .vcf <- subset(.vcf, nchar(from) == 1 &</pre>
                   nchar(to) == 1 & !is.na(qual)) # no indels
  .infosplit <- str_split_fixed(.vcf$info, regex("[=;]"), 15)[,c(2, 4, 12, 14)]</pre>
  .vcf$coverage <- as.numeric(.infosplit[,1])</pre>
```

```
.vcf$frequency <- as.numeric(.infosplit[,2])</pre>
  .vcf$bb.a <- as.numeric(.infosplit[,3])</pre>
  .vcf$bb.b <- as.numeric(.infosplit[,4])</pre>
  .vcf$info <- NULL</pre>
  .vcf$count <- as.integer(round(.vcf$coverage * .vcf$frequency))</pre>
  .vcf
}
head(read_vcf("p126.h4_2_ballast_m1.vcf"))
## Warning in read_vcf("p126.h4_2_ballast_m1.vcf"): NAs introduced by coercion
##
     chromosome position from to qual coverage
                                                      frequency
## 1
           chr2 212295704
                                              2265 0.0004415011 0.9957058
                              C A
                                       8
## 2
           chr2 212295705
                              C A
                                       8
                                              2265 0.0004415011 0.9957058
## 3
                              A G
                                       7
                                              2265 0.0004415011 1.0295689
           chr2 212295713
## 4
           chr2 212295718
                              G A
                                       6
                                              2265 0.0004415011 2.1314828
## 5
           chr2 212295725
                              C A
                                      8
                                              2264 0.0004416961 0.9957058
## 6
           chr2 212295732
                              C T
                                      24
                                              2264 0.0013250883 2.1314828
         bb.b count
## 1 21862.23
## 2 21862.23
                   1
## 3 16394.12
                   1
## 4 24344.82
                   1
## 5 21862.23
                   1
## 6 24344.82
Read samples with HD734 standard DNA and control human DNA, append metadata
df <- data.frame()</pre>
read_vcf_with_metadata <- function(file_name, primer_set, replica, ratio, type) {</pre>
  .vcf <- read vcf(file name)</pre>
  .vcf <- merge(.vcf, df.vmeta, all.x = type != "standard", all.y = F)</pre>
  .vcf$known.frequency <- .vcf$known.frequency * ratio</pre>
  .vcf$known.frequency[is.na(.vcf$known.frequency)] <- 0</pre>
  .vcf$primer_set <- primer_set</pre>
  .vcf$replica <- primer_set</pre>
  .vcf$type <- type</pre>
  .vcf <- subset(.vcf, frequency < 0.4 & count > 0) # remove alleles in control
  .vcf
for (i in 1:nrow(df.smeta)) {
  df <- with(df.smeta, rbind(df,</pre>
                               read_vcf_with_metadata(paste(prefix[i], "vcf", sep="."),
                                                       primer_set[i],
                                                       replica[i],
                                                       ratio[i],
                                                       type[i])))
```

```
## Warning in read_vcf(file_name): NAs introduced by coercion
## Warning in read_vcf(file_name): NAs introduced by coercion
## Warning in read_vcf(file_name): NAs introduced by coercion
## Warning in read vcf(file name): NAs introduced by coercion
## Warning in read_vcf(file_name): NAs introduced by coercion
Group observed variants into tiers
df$tier <- cut(df$known.frequency, c(-1, 0, 0.009, 0.02, 1))
levels(df$tier) <- c("error", "0.1%", "1%", "5%")</pre>
summary(df$tier)
         0.1%
                  1%
                        5%
## error
    2468
                  46
                        12
Number of detected variants and sampling
df.v.summary <- df.vmeta</pre>
df.v.summary$known.frequency <- df.v.summary$known.frequency * 0.1 # 1:9 diluted samples
df.v.summary$ratio <- "1:10"</pre>
df.v.summary.2 <- df.vmeta</pre>
df.v.summary.2$ratio <- "1:1"</pre>
```

df.v.summary\$tier <- cut(df.v.summary\$known.frequency, c(-1, 0, 0.009, 0.02, 1))

df.v.summary <- rbind(df.v.summary, df.v.summary.2)</pre>

```
levels(df.v.summary$tier) <- c("error", "0.1%", "1%", "5%")
print(df.v.summary)</pre>
```

##		chromosome	position	from	
##	1	chr7	140453136	A	T
##	2	chr7	140453136	A	C
##	3	chr7	140453137	C	T
##	4	chr7	55241707	G	A
##	5	chr7	55259515	T	G
##	6	chr7	55259524	T	A
##	7	chr7	55249071	C	T
##	8	chr19	3118942	A	T
##	9	chr2	209113112	C	T
##	10	chr2	209113113	G	A
##	11	chr12	25378562	C	T
##	12	chr12	25398284	C	G
##	13	chr12	25398285	C	A
##	14	chr12	25398284	C	T
##	15	chr12	25398285	C	G
##	16	chr12	25398285	C	T
##	17	chr12	25398284	C	A
##	18	chr12	25398281	C	T
##	19	chr12	25380275	T	G
##	20	chr12	25380276	T	A
##	21	chr1	115256528	T	A
##	22	chr1	115256530	G	T
##	23	chr1	115256529	T	A
##	24	chr1	115256529	T	C
##	25	chr3	178936082	G	A
##	26	chr3	178936091	G	A
##	27	chr3	178952085	A	G
##	28	chr7			G
##	29	chr7	140453136	A	T
##	30	chr7	140453136	A	C
##	31	chr7	140453137	C	T
##	32 33	chr7	55241707	G	A
##		chr7	55259515 55259524	T	G
##	34	chr7	55249071	T	A
##	35	chr7 chr19		C	T T
##	36		3118942	A	
	37		209113112	C	T
##		chr2		G	A
##	39	chr12	25378562	C	T
	40	chr12	25398284	C	G
	41	chr12	25398285	C	A
	42	chr12 chr12	25398284	C	T
	43		25398285	C	G
	44	chr12	25398285	C	T
	45 46	chr12	25398284	C	A
		chr12	25398281		T
	47	chr12	25380275	T	G ^
	48	chr12	25380276	T	A
##	49	chr1	115256528	Т	A

```
## 50
             chr1 115256530
## 51
                                             Т
             chr1 115256529
                                                Α
## 52
             chr1 115256529
                                             Τ
             chr3 178936082
                                             G
## 53
                                                Α
## 54
             chr3 178936091
                                             G
                                                Α
                                                G
## 55
             chr3 178952085
                                             Α
## 56
             chr7
                   55242465 GGGAATTAAGAGAAGC
                                                G
##
                                       id known.frequency
                                                            reference ratio tier
## 1
                  BRAF_V600E, BRAF_V600K
                                                   0.0080
                                                             BRAF_E15
                                                                       1:10 0.1%
##
                  BRAF_V600G, BRAF_V600R
                                                   0.0010
                                                             BRAF_E15
                                                                       1:10 0.1%
      BRAF_V600M, BRAF_V600K, BRAF_V600R
                                                   0.0010
                                                             BRAF_E15
                                                                        1:10 0.1%
                                                             EGFR_E18
## 4
                              EGFR_G719S
                                                   0.0167
                                                                        1:10
## 5
                              EGFR_L858R
                                                   0.0010
                                                             EGFR_E21
                                                                       1:10 0.1%
                                                                        1:10 0.1%
## 6
                              EGFR_L861Q
                                                   0.0010
                                                             EGFR_E21
## 7
                              EGFR_T790M
                                                   0.0010
                                                             EGFR_E20
                                                                        1:10 0.1%
## 8
                             GNA11_Q209L
                                                   0.0013
                                                             GNA11_E5
                                                                        1:10 0.1%
## 9
                                                   0.0013
                                                              IDH1_E4
                              IDH1_R132C
                                                                        1:10 0.1%
## 10
                              IDH1 R132H
                                                   0.0013
                                                              IDH1 E4
                                                                       1:10 0.1%
## 11
                              KRAS_A146T
                                                   0.0013
                                                              KRAS_E4
                                                                       1:10 0.1%
## 12
                               KRAS_G12A
                                                   0.0013
                                                              KRAS E2
                                                                       1:10 0.1%
## 13
                               KRAS_G12C
                                                   0.0013
                                                              KRAS_E2
                                                                       1:10 0.1%
                               KRAS_G12D
                                                   0.0013
                                                              KRAS E2
## 14
                                                                       1:10 0.1%
                               KRAS_G12R
                                                              KRAS_E2
## 15
                                                   0.0013
                                                                       1:10 0.1%
## 16
                               KRAS G12S
                                                   0.0013
                                                              KRAS E2
                                                                       1:10 0.1%
## 17
                               KRAS_G12V
                                                   0.0013
                                                              KRAS E2
                                                                       1:10 0.1%
## 18
                               KRAS_G13D
                                                   0.0250
                                                              KRAS_E2
                                                                       1:10
                                                                               5%
## 19
                               KRAS_Q61H
                                                   0.0013
                                                              KRAS_E3
                                                                        1:10 0.1%
## 20
                               KRAS_Q61L
                                                   0.0013
                                                              KRAS_E3
                                                                       1:10 0.1%
## 21
                                                              NRAS_E3
                               NRAS_Q61H
                                                   0.0013
                                                                       1:10 0.1%
## 22
                               NRAS_Q61K
                                                   0.0013
                                                              NRAS_E3
                                                                       1:10 0.1%
## 23
                               NRAS_Q61L
                                                   0.0013
                                                              NRAS_E3
                                                                        1:10 0.1%
## 24
                               NRAS_Q61R
                                                   0.0013
                                                              NRAS_E3
                                                                        1:10 0.1%
## 25
                           PIK3CA_E542K
                                                   0.0013
                                                            PIK3CA_E9
                                                                        1:10 0.1%
## 26
                           PIK3CA_E545K
                                                   0.0013
                                                            PIK3CA_E9
                                                                        1:10 0.1%
## 27
                          PIK3CA H1047R
                                                   0.0300 PIK3CA E20
                                                                        1:10
                                                                               5%
## 28
                        EGFR_dE746-A750
                                                   0.0010
                                                             EGFR_E19
                                                                        1:10 0.1%
## 29
                  BRAF V600E, BRAF V600K
                                                   0.0800
                                                             BRAF E15
                                                                         1:1
                  BRAF_V600G, BRAF_V600R
                                                   0.0100
                                                             BRAF_E15
                                                                         1:1
                                                                                1%
## 30
      BRAF_V600M, BRAF_V600K, BRAF_V600R
                                                             BRAF_E15
                                                                                1%
## 31
                                                   0.0100
                                                                         1:1
## 32
                                                                               5%
                              EGFR_G719S
                                                   0.1670
                                                             EGFR_E18
                                                                         1:1
## 33
                              EGFR L858R
                                                   0.0100
                                                             EGFR E21
                                                                         1:1
                                                                               1%
## 34
                                                                               1%
                              EGFR L861Q
                                                   0.0100
                                                             EGFR E21
                                                                         1:1
## 35
                              EGFR T790M
                                                   0.0100
                                                             EGFR E20
                                                                         1:1
                                                                               1%
## 36
                                                                               1%
                             GNA11_Q209L
                                                   0.0130
                                                             GNA11_E5
                                                                         1:1
## 37
                              IDH1_R132C
                                                   0.0130
                                                              IDH1_E4
                                                                         1:1
                                                                               1%
## 38
                                                              IDH1_E4
                                                                                1%
                              IDH1_R132H
                                                   0.0130
                                                                         1:1
## 39
                              KRAS_A146T
                                                   0.0130
                                                              KRAS_E4
                                                                         1:1
                                                                               1%
## 40
                               KRAS_G12A
                                                   0.0130
                                                              KRAS_E2
                                                                         1:1
## 41
                               KRAS_G12C
                                                   0.0130
                                                              KRAS_E2
                                                                         1:1
                                                                               1%
                                                                         1:1
## 42
                               KRAS_G12D
                                                   0.0130
                                                              KRAS_E2
                                                                               1%
## 43
                                                              KRAS_E2
                               KRAS_G12R
                                                   0.0130
                                                                         1:1
                                                                               1%
## 44
                               KRAS_G12S
                                                   0.0130
                                                              KRAS_E2
                                                                         1:1
                                                                               1%
## 45
                               KRAS_G12V
                                                   0.0130
                                                              KRAS_E2
                                                                         1:1
                                                                               1%
## 46
                               KRAS G13D
                                                   0.2500
                                                              KRAS E2
                                                                               5%
```

```
## 47
                             KRAS Q61H
                                                 0.0130
                                                           KRAS E3
                                                                     1:1
                                                                           1%
## 48
                             KRAS_Q61L
                                                           KRAS E3
                                                                     1:1
                                                                            1%
                                                 0.0130
## 49
                             NRAS Q61H
                                                 0.0130
                                                           NRAS E3
                                                                     1:1
                                                                            1%
## 50
                             NRAS_Q61K
                                                 0.0130
                                                           NRAS_E3
                                                                     1:1
                                                                            1%
## 51
                             NRAS Q61L
                                                 0.0130
                                                           NRAS E3
                                                                     1:1
                                                                            1%
## 52
                             NRAS Q61R
                                                 0.0130
                                                           NRAS E3
                                                                     1:1
                                                                           1%
## 53
                          PIK3CA E542K
                                                 0.0130 PIK3CA E9
                                                                     1:1
                                                                           1%
                                                 0.0130 PIK3CA E9
## 54
                          PIK3CA E545K
                                                                     1:1
                                                                           1%
## 55
                         PIK3CA H1047R
                                                 0.3000 PIK3CA E20
                                                                     1:1
                                                                           5%
                                                          EGFR_E19
## 56
                       EGFR_dE746-A750
                                                 0.0100
                                                                     1:1
                                                                           1%
df.v.summary <- df.v.summary %>%
  group_by(tier) %>%
  summarise(total = n())
df.v.summary$total <- 2 * df.v.summary$total # we have 2 independent experiments
df.v.summary <- merge(df.v.summary,</pre>
                      df %>% group_by(tier) %>%
                        summarise(found = n(),
                                  detect.prob.m = mean(1-exp(-frequency * coverage)),
                                  detect.prob.sd = sd(1-exp(-frequency * coverage))))
df.v.summary$detection.prob.lb <- with(df.v.summary,</pre>
                                        detect.prob.m - 1.96*detect.prob.sd - 5e-3)
df.v.summary$detection.prob.ub <- with(df.v.summary,</pre>
                                       detect.prob.m + 1.96*detect.prob.sd)
df.v.summary$detection.prob.ub <- ifelse(df.v.summary$detection.prob.ub > 1,
                                          df.v.summary$detection.prob.ub )
library(scales)
fig3 <- ggplot(df.v.summary, aes(x=tier)) +</pre>
  geom_bar(stat="identity", aes(y=detection.prob.ub), fill = "grey70", width=0.5)+
  geom_bar(stat="identity", aes(y=detection.prob.lb), fill = "white")+
  geom_point(aes(y = found/total), size = 5, color="red", shape="+") +
  geom_text(aes(y = found/total - 0.05, label = paste(found, total, sep="/"))) +
  scale_y_continuous("Detection rate", labels=percent) +
  xlab("Variant tier") +
  theme bw() +
  theme(panel.grid.major.x = element_blank(),
        panel.grid.major.y = element_blank())
save(file = "../figures/fig3.Rda", fig3)
fig3
Quality score and error P-values
df.e <- df %>%
 filter(type == "blank") %>%
  mutate(true.p.value = -10 * log10(1 - (rank(frequency) - 0.5) / n()))
```

fig4 <- ggplot(df.e, aes(true.p.value, qual)) +</pre>




Figure 1: Number of detected variant for each variant frequency tier across two independent experiments with the reference standard. Shaded areas show the 95% confidence intervals for expected fraction of recovered variants, i.e. binomial proportion confidence intervals built using known variant frequency and template coverage.

```
geom_abline(slope = 1, intercept = 0, linetype = "dashed") +
geom_smooth(color="red", fill="grey", alpha=0.5) +
scale_x_continuous("-10 log10 P-value", limits = c(0, 40)) +
scale_y_continuous("Q score", limits = c(0, 40)) +
theme_bw()
save(file = "../figures/fig4.Rda", fig4)
```

fig4

`geom_smooth()` using method = 'gam'

Warning: Removed 1 rows containing non-finite values (stat_smooth).




Figure 2: MAGERI Q score and the empirical P-values of erroneous variants detected in control donor DNA.

```
nrow(df.e)
## [1] 2468
cor(df.e$qual, df.e$true.p.value, method="pearson")
## [1] 0.8311097
Frequency of variants and errors
fig5 <- ggplot(df, aes(x=frequency, color = tier)) +
    geom_density(aes(fill=tier), alpha=0.7) +
    ylab("Density") +</pre>
```

```
scale_x_log10("Frequency") +
scale_color_brewer("Variant tier", palette = "Set1") +
scale_fill_brewer("Variant tier", palette = "Set1") +
theme_bw()
save(file = "../figures/fig5.Rda", fig5)
```

fig5




Figure 3: Frequency distribution of known Tru-Q 7 variants coming from each frequency tier and errors in the control donor DNA.

```
df.known <- subset(df, known.frequency > 0)
nrow(df.known)

## [1] 101
mean((df.known$frequency - df.known$known.frequency) / df.known$known.frequency)

## [1] 0.07155274

cor(df.known$frequency, df.known$known.frequency, method="spearman")

## [1] 0.830348

Variant quality scores

fig6 <- ggplot(df,aes(tier, qual, fill=tier)) +
    geom_hline(yintercept = 13.0103, linetype = "dotted") +</pre>
```

```
geom_hline(yintercept = 20.0103, linetype ="dashed") +
geom_quasirandom(varwidth = T, shape=21, color="grey10") +
annotate("text", label = "p < 0.05", x=4.3,y=13.0103+2) +
annotate("text", label = "p < 0.01", x=4.3,y=20.0103+4) +
scale_fill_brewer("", palette = "Set1", guide=F) +
xlab("Variant tier") + ylab("Q score") +
theme_bw()
save(file = "../figures/fig6.Rda", fig6)</pre>
```

fig6




Figure 4: Comparison of Q score distribution of erroneous variants and variants of each frequency tier. Dotted and dashed lines show P < 0.05 and P < 0.01 thresholds respectively.

ROC curve for rare variants

```
df.1 <- subset(df, tier %in% c("error", "0.1%"))
df.1$type <- ifelse(df.1$tier == "error", 0, 1)

make_roc <- function(rocobj, type) {
   .df.roc <- data.frame(spec = seq(0, 1, 0.01))
   sens.ci <- ci.se(rocobj, specificities = .df.roc$spec)

   .df.roc$sens.lo <- sens.ci[,1]
   .df.roc$sens.me <- sens.ci[,2]
   .df.roc$sens.hi <- sens.ci[,3]</pre>
```

```
.df.roc$type <- type</pre>
  .df.roc
}
rocobj <- roc(type ~ qual, df.1, ci=T)</pre>
print(rocobj)
##
## Call:
## roc.formula(formula = type ~ qual, data = df.1, ci = T)
## Data: qual in 2468 controls (type 0) < 43 cases (type 1).
## Area under the curve: 0.9277
## 95% CI: 0.8727-0.9827 (DeLong)
df.roc <- make_roc(rocobj, "Q score")</pre>
rocobj <- roc(type ~ frequency, df.1, ci=T)</pre>
print(rocobj)
##
## Call:
## roc.formula(formula = type ~ frequency, data = df.1, ci = T)
## Data: frequency in 2468 controls (type 0) < 43 cases (type 1).
## Area under the curve: 0.8637
## 95% CI: 0.7839-0.9435 (DeLong)
df.roc <- rbind(df.roc, make_roc(rocobj, "Frequency"))</pre>
fig7 <- ggplot(df.roc, aes(x=spec)) +</pre>
  geom_ribbon(aes(ymin=sens.lo, ymax=sens.hi, group=type), fill="grey", alpha=0.5) +
  geom_abline(slope = 1, intercept = 1, linetype = "dashed") +
  geom_line(aes(y=sens.me, color = type)) +
  scale_x_reverse("Specificity") +
  scale_y_continuous("Sensitivity", limits=c(0,1)) +
  scale_color_brewer("Threshold", palette = "Set1") +
 theme bw()
save(file = "../figures/fig7.Rda", fig7)
fig7
```




Figure 5: Receiver operation characteristic (ROC) curve comparing the sensitivity and specificity of MAGERI Q scores (blue line) and frequency-based thresholding (red line) in the task of classification of errors and 0.1% tier variants.